

D95+ Surgical Respirator Mask

The D95+ is a single-use surgical respirator mask. You must be Fit Tested prior to using the D95+ surgical respirator mask.

Please contact your local Network Staff Health Service, Worker Health Team or SAAS fit test coordinator or for Regional Networks, your Infection Control Nurse to organise a Fit Test appointment.

Donning and Doffing Instructions

Donning – correctly fitting the D95+ respirator mask



1. Prior to donning, practice hand hygiene.
2. Ensure you select the correct size D95+ mask that you have been fitted for (small, regular or large).
3. Move the Fitclip™ to the band on your preferred side for donning*
4. Open the mask, flattening the nose piece.
5. Lean chin into the mask and lift bands toward the crown of the head, securing onto the Fitclip™.
6. Adjust mask bands for a secure and comfortable fit.

Now, you must perform a fit check



Fit Check – Must be performed each time a D95+ respirator mask is donned.

1. Mould the nose wire securely. Pinch the nose piece (angled over the nostril) and drop the jaw gently to take the natural shape of the face.
2. Look left, right, up and down. Ensure there is NO movement in the mask at the nose piece.
3. Placing hands above the mask, deeply inhale and exhale feeling for escape of air.

There should be no gaps between the face and mask, providing an adequate seal. Once the mask is applied it should not be touched or reapplied.

If you have not achieved a successful fit as instructed above you MUST seek advice from Staff Health or Worker Health Nurse. Regional LHN staff will need to contact Infection Control. An incorrectly fitted respirator mask may not provide you with the intended level of protection.

Doffing – safely removing the mask

1. Practice hand hygiene
2. Without touching the mask, follow bands towards FitClip™ and gently unhook, keeping tension on the bands.
3. Lean forward maintaining control of the band tension, whilst slowly removing the mask away from the face.
4. Ensure the mask does not flick forward or back onto the face.
5. Discard safely.
6. Practice hand hygiene.



*FitClip™ design may vary.